

Item No. 92230

one step disinfectant cleaning wipes

CLEANS, SANITIZES, DISINFECTS & DEODORIZES
***VIRUCIDAL / BACTERICIDAL / PSEUDOMONICIDAL**

- Kills 99.9% of Bacteria in just 15 Seconds**
- Kills Antibiotic Resistant Bacteria: CA-MRSA, MRSA, VISA, VRE and NDM-1
- Kills Cold¹ & Flu Viruses^{†††}
- Disinfects in just 5 minutes²
- Orange Scented
- The smart way to clean and sanitize and disinfect



ACTIVE INGREDIENTS

Alkyl (60% C14, 30% C16, 5% C12, 5% C18)	
dimethyl benzyl ammonium chlorides	0.14%
Alkyl (68% C12, 32% C14) dimethyl ethylbenzyl ammonium chlorides	0.14%
INERT INGREDIENTS	99.72%
TOTAL	100.00%

**KEEP OUT OF REACH OF CHILDREN
CAUTION**

See back panels for additional precautionary statements.



EPA Reg. No. 1839-190-70925 EPA Est. No. 70925-CA-001
NET CONTENTS: 2 lbs. 0.45 ozs., 130 Pre-moistened Wipes (6"x8")
MADE IN USA

One Step Disinfectant Cleaning Wipes are intended for use in:

- Hospitals and Healthcare settings: ambulances, anesthesia, CAT labs, CCU, Central Supply, Clinics, Dental Offices, Dialysis Clinics, Doctor's offices, Emergency Rooms, E.R., Home health care, Hospices, ICU, Laboratory, Laundry rooms, Neonatal Intensive Care Units (NICU), Newborn nursery, Nursing homes, Operating rooms, Pediatrics, Pediatric Intensive Care Units (PICU), Physical therapy, Ophthalmic offices, Radiology, Recovery Room, Rehabilitation, Respiratory therapy, Surgery, Transport vehicles, X-Ray
- Veterinary clinics, animal hospitals, animal life science laboratories, kennels, breeding establishments, grooming establishments, pet animal quarters, zoos, pet shops, tack shops
- Beauty salons, barber shops, correctional facilities, funeral homes, mortuaries
- Schools, health clubs, gyms, athletic facilities, airports, transportation terminals
- Grocery stores, Shopping Centers, Delicatessens, Catering facilities, Banquet halls
- Nursery, Daycare centers
- Bathrooms
- Kitchens[#], Restaurants[#]

Use One Step Disinfectant Cleaning Wipes to clean and disinfect non-porous personal protective safety equipment, protective headgear, hard hats, half mask respirators, full face breathing apparatus, gas masks, goggles, spectacles, face shields, silicone rubber or PVC hearing protectors and vinyl covered ear muffs. Rinse all equipment that comes in prolonged contact with skin before reuse with clean warm water (about 120°F), and allow to air dry. (Caution: Cleaning at the 120°F temperature will avoid overheating and distortion of the personal safety equipment that would necessitate replacement.)

One Step Disinfectant Cleaning Wipes will clean and disinfect hard, non-porous surfaces of: cabinets, chairs, computers, countertops, desks, doorknobs, fax machines, headsets, keyboards, railings, seats, sinks, tables, telephones, trashcans, appliance exteriors, faucets, bathtubs, shower stalls, stove tops, grocery carts, changing tables, cribs, diaper pails, gym equipment, exercise machines, wrestling mats, gymnastic equipment.

DIRECTIONS FOR USE

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. This product is not for use on critical and semi-critical medical device surfaces.

CLEANING: Use one One Step Disinfectant Cleaning Wipes to completely remove all gross dirt and debris from surfaces to be treated.

DISINFECTING: To disinfect hard, non-porous surfaces, thoroughly wet the surface to be treated. Treated surface must remain visibly wet for disinfecting. Use enough wipes for treated surface to remain visibly wet for 5 minutes. For heavily soiled surfaces, clean before following disinfecting instructions. Allow surfaces to air dry.

***VIRUCIDAL ACTIVITY:** This product is an effective virucide on hard, non-porous surfaces against Influenza A2/Japan, Herpes Simplex Virus 1 (ATCC VR-260), Human Immunodeficiency Virus, Type 1 (HIV-1) and Herpes Simplex Virus 2 (ATCC VR-734), Duck Hepatitis B Virus (DHBV), Human Hepatitis B Virus (HBV), Bovine Viral Diarrhea Virus (BVDV) and Human Hepatitis C Virus (HCV) when the treated surface is allowed to remain wet for 5 minutes. Allow surface to air dry.

¹This product is an effective virucide on hard, non-porous surfaces against Human Coronavirus (ATCC VR-740) when the treated surface is allowed to remain wet for 10 minutes. Allow surfaces to air dry.

KILLS HIV-1 ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS in health care settings (Hospitals, Nursing Homes) or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1, Human Hepatitis B Virus (HBV) and Human Hepatitis C Virus (HCV).

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1, HBV and HCV OF SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS.

PERSONAL PROTECTION: When handling items soiled with blood or body fluids, use disposable latex gloves, gowns, masks and eye coverings.

CLEANING PROCEDURES: Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of this product.

CONTACT TIME: Allow surface to remain wet for 5 minutes.

DISPOSAL OF INFECTIOUS MATERIALS: Blood and other body fluids should be autoclaved and disposed of according to local regulations for infectious waste disposal.

FUNGICIDAL ACTIVITY: This product is an effective fungicide against Trichophyton mentagrophytes (ATCC 9533) (causative agent of Athlete's Foot Fungus) (a cause of ringworm) when used as directed on hard surfaces found in bathrooms, shower stalls, locker rooms, exercise facilities or other clean, hard, non-porous, surfaces commonly

contacted by bare feet. The surface must remain visibly wet for 3 minutes. Allow surfaces to air dry.

EFFICACY TESTS HAVE DEMONSTRATED THAT THIS PRODUCT IS AN EFFECTIVE BACTERICIDE, VIRUCIDE AND FUNGICIDE IN THE PRESENCE OF ORGANIC SOIL (5% BLOOD SERUM).

SANITIZING/DEODORIZING: To sanitize/deodorize hard, non-porous surfaces; wipe surface, using enough wipes to allow treated surface to remain wet for 15 seconds for Staphylococcus aureus (Staph) (ATCC 6538), Klebsiella pneumoniae (ATCC 4352), Listeria monocytogenes (ATCC 35152), Escherichia coli (E. coli) O157:H7 (ATCC 43895) and 60 seconds for Campylobacter jejuni (ATCC 29428). Let surface dry. For highly soiled surfaces, clean excess dirt first.

For surfaces that may come in contact with food, a potable water rinse is required. For plastic and painted surfaces, spot test on an inconspicuous area before use. This product will not harm most surfaces including acrylic, fiberglass and vinyl. **Do not use on unpainted wood.** Not for personal cleansing. This is **not** a baby wipe! Do not flush in toilet.

TO OPEN PACKAGE: Flip open dispensing cap and remove lid by pulling upward. Thread center wipe corner through slits in lid center. DO NOT PUSH FINGER THROUGH OPENING. Replace lid. Dispense wipes by pulling out at an angle. Snap flip lid down when finished, to retain moisture. Dispose of used wipe in trash. Do not flush.

STORAGE AND DISPOSAL

Do not contaminate water, food or feed by storage or disposal.
CONTAINER STORAGE/DISPOSAL: Store in an area inaccessible to children. Wipes: Dispose of wipe in trash after use. Do not flush. Canister: Tightly close lid between uses to retain moisture. Nonrefillable container. Do not reuse or refill this canister. Offer empty canister for recycling if available or discard in trash.

PRECAUTIONARY STATEMENTS

Hazards to Humans and Domestic Animals

CAUTION: Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling.

FIRST AID

Have the product container or label with you when calling a poison control center or doctor, or going for treatment.

IF IN EYES: Hold eyes open and rinse slowly and gently with water for 15–20 minutes. Remove contact lenses if present, after the first 5 minutes, then continue rinsing. Call a poison control center or doctor for treatment advice.

**Staphylococcus aureus (Staph), Klebsiella pneumoniae, Listeria monocytogenes, Escherichia coli (E. coli) O157:H7
 †††Influenza A Virus, Pandemic 2009 H1N1 Influenza A Virus, Human Coronavirus¹
²This product is bactericidal in the presence of organic soil (5% blood serum) against Escherichia coli (E. coli) O157:H7 (ATCC 43895), Salmonella enterica (ATCC 10708), Pseudomonas aeruginosa (ATCC 15442), Staphylococcus aureus (Staph) (ATCC 6538), Streptococcus pyogenes (Strep) (ATCC 19615), Community Associated Methicillin resistant Staphylococcus aureus Genotypes 300 and 400 (CA-MRSA), Methicillin resistant Staphylococcus aureus (MRSA), Vancomycin resistant Enterococcus faecium (VRE), Vancomycin intermediate resistant Staphylococcus aureus (VISA) and Klebsiella pneumoniae New Delhi Metallo-Beta Lactamase (NDM-1) Carbapenem Resistant (CDC Clinical Isolate No. 10002) at 5 minutes contact time.